

# 一次性使用采样拭子尺寸详情

简石生物  
JIANSHI BIOTECH

## 咽拭子



拭子总长度  $150 \pm 1 \text{ mm}$   
折断点长度  $30 \pm 1 \text{ mm}$



拭子总长度  $150 \pm 1 \text{ mm}$   
折断点长度  $76 \pm 1 \text{ mm}$



拭子总长度  $150 \pm 1 \text{ mm}$   
折断点长度  $20 \pm 1 \text{ mm}$



拭子总长度  $150 \pm 1 \text{ mm}$   
折断点长度  $0 \text{ mm}$

## 鼻拭子



拭子总长度  $150 \pm 1 \text{ mm}$   
折断点长度  $80 \pm 1 \text{ mm}$



拭子总长度  $90 \pm 1 \text{ mm}$   
折断点长度  $50 \pm 1 \text{ mm}$

## 宫颈拭子



拭子总长度  $170 \pm 1 \text{ mm}$   
折断点长度  $50 \pm 1 \text{ mm}$



拭子总长度  $170 \pm 1 \text{ mm}$   
折断点长度  $90 \pm 1 \text{ mm}$