

一次性使用采样拭子尺寸详情

简石生物
JIANSHI BIOTECH

咽拭子



拭子总长度 150 ± 1 mm
折断点长度 30 ± 1 mm



拭子总长度 150 ± 1 mm
折断点长度 76 ± 1 mm



拭子总长度 150 ± 1 mm
折断点长度 20 ± 1 mm



拭子总长度 150 ± 1 mm
折断点长度 0 mm

鼻拭子



拭子总长度 150 ± 1 mm
折断点长度 80 ± 1 mm



拭子总长度 90 ± 1 mm
折断点长度 50 ± 1 mm

宫颈拭子



拭子总长度 150 ± 1 mm
折断点长度 50 ± 1 mm



拭子总长度 150 ± 1 mm
折断点长度 90 ± 1 mm